

MotoRater

Ð

MotoRater

Innovative Kinematic Movement & Gait Analysis for Mice & Rats

1

www.TSE-Systems.com | Info@TSE-Systems.com

MotoRater

MotoRater "MotionRater" is an innovative automated tool for standardized quantitative and objective evaluation of fine motor function and kinematic analysis in small laboratory animals.

Movement analysis is not confined to footprints; it also analyzes kinematic parameters for different motion endpoints. The MotoRater observes the animal from three sides simultaneously and provides explicit readouts for relevant body parts, i.e., limbs, joints, tail, etc. making this a comparative translational apparatus for disease study to humans. The modular setup allows evaluation from all kinematic movements during four-movement modalities (overground locomotion, skilled ladder walking, wading, and swimming) providing a highly differentiated analysis of the animal's motor performance.

The option to test animals in the wading and swimming module provides the unique opportunity to evaluate motion patterns of severely impaired animals, which would not support their body weight on solid around.



Benefits

Automated kinematic movement and gait analysis

A single system tests four different movement modalities

Virtually unlimited Simultaneous observation and number of analysis of three parameters animal sides

Suitable for many (inclunding severely impaired) animal models

TSE Motion powered by Simi Reality Motion Systems is software for video acquisition and data analysis.

To make it simpler, the NEW MotoRater can be operated directly through one software for video acquisition and data analysis providing video through a stationary high-speed "fish-eye objective" camera and analysis through a marker-based kinematic gait analysis for maximal variability and scientific analysis.

TSE Motion: Video Acquisition

- The animal is placed into the activity compartment and freely moves through the long activity runway (either dry, on a ladder, or containing water) towards a home cage or goal box (serving as a motivator); the long runway favors natural gait and allows recording of 3-4 gait cycles.
- Video properties can be adjusted, permitting optimal tracking performance for various acquisition conditions and animal color.



TSE Motion Analysis software tracks any type and number of specified body points (automatized tracking) during four different movement modalities:



- Kinematically meaningful connections between two or more body points can be defined for calculating joint angles, distances, velocity, and acceleration, etc.
- Body points can be tracked in all frames or only a selection of relevant frames.
- Automated or manual point-and-click tracking is also possible. (relevant if only one or a few selected frames are of interest).

More than 100 kinematic gait parameters

- Exact & standardized analysis protocols can be created, saved, and applied to tracking data maximizing data reliability across animals and minimizing group size
- Objective and unbiased analysis
- Uniquely comprehensive calculation options (including distances, angles, mean, maximum and minimum values, etc.) for:

General spatiotemporal patterns: Inter-limb coordination: Swing phase and paw trajectory: Body posture and joint angles: Gait variability, deviations:

retraction, protraction

- Data smoothing and filtering options
- Elaborate options for graphical data representation including stick view figures, phase model diagrams, statistic plots, and more
- Dynamic video and diagram animation
- Programming of analysis templates; standard templates (e.g. gait analysis) are available on request

Research Applications

Animals are objectively and reliably screened for deficits in a high-throughput fashion.

Relevant disease models include:

- Ischemia / Stroke
- Huntington's disease
- Multiple sclerosis
- Amyotrophic lateral sclerosis
- Duchenne muscular dystrophy
- Parkinson's disease
- Batten disease

• Pain / Arthritis



stride distance, duration, speed, stance time, swing time

- diagonal cadence, left/right alternation rhythm
- swing speed, smoothness, trajectory shape, abnormality
- hip height, tail & nose height, limb function

 Spinal cord injury Spasticity / Ataxia

 Locomotor recovery / treatment / exercise efficacy

Other neurological and rare diseases



TSE MotoRater can be integrated with other cutting-edge technologies

As a synchronization hub, the trigger box allows advanced users to integrate external systems into the MotoRater. For example, EMG recordings can be triggered and observed directly in the TSE Motion software. Additionally, the wireless EMG option can be obtained through TSE Stellar Telemetry and wireless optogenetics.



Selected Publications

The MotoRater from TSE Systems is world-renowned among leading academic institutes, research centers, contract research organizations, and pharmaceutical companies. Dozens of scientific papers based their research on the MotoRater have been published in peer-reviewed journals like Nature, Cell, and others. You can find below a selection of scientific publications.

Yoshizaki S, Kijima K, Hara M, Saito T, Tamaru T, Tanaka M, Konno D-J, Nakashima, Okada S. **Tranexamic acid reduces heme** cytotoxicity via the TLR4/TNF axis and ameliorates functional recovery after spinal cord injury. J Neuroinflammation. 2019; 16: 1–15.

Krupa P, Svobodova B, Dubisova J, Kubinova S, Jendelova P, Machova Urdzikova L. Nano-formulated curcumin (LipodisqTM) modulates the local inflammatory response, reduces glial scar, and preserves the white matter after spinal cord injury in rats. Neuropharmacology. 2019; 155: 54–64.

Chen B, Li Y, Yu B, Zhang Z, Brommer B, Williams PR, Liu Y, Hegarty SV, Zhou S, Zhu J, Guo H, Lu Y, Zhang Y, Gu X, He Z. Reactivation of Dormant Relay Pathways in Injured Spinal Cord by KCC2 Manipulations. Cell. 2018; 174: 521-535.

Preisig DF, Kulic L, Krüger M, Wirth F, McAfoose J, Späni C, Gantenbein P, Derungs R, Nitsch RM, Welt T. **High-speed video gait** analysis reveals early and characteristic locomotor phenotypes in mouse models of neurodegenerative movement disorders. Behav Brain Res. 2016; 311: 340–353.

Zörner B, Filli L, Starkey ML, Gonzenbach R, Kasper H, Röthlisberger M, Bolliger M, Schwab ME. **Profiling locomotor recovery:** Comprehensive quantification of impairments after CNS damage in rodents. Nature Methods. 2010; 7 (9): 701-708.

<u>Novel Principle Component Analysis to Assess Gait in Chronically Exercised vs. Unexercised MDX Mice Shows Both Exacerbation</u> and Amelioration of the Underlying Phenotype (POSTER)

Fine Motor Kinematic Analysis

The information in this document is in good faith, and while every care has been taken in preparing these documents, TSE Systems makes no representations and gives no warranties of whatever nature in respect of these documents, including but not limited to the accuracy or completeness of any information, facts and/or opinions contained therein.

Copyright © 2021 TSE Systems - All rights reserved.