



PhenoMaster NG

One platform. Complete phenotyping.



The Next-Generation Home-Cage Phenotyping for Rodents

As the most advanced home-cage phenotyping platform, PhenoMaster enables researchers to integrate metabolic, behavioral, and physiological measurements into one synchronized system – continuously, automatically, and non-invasively.

Its modular architecture covers every aspect of in vivo research – from indirect calorimetry and activity monitoring to microbiome interactions and feeding behavior analysis – generating precise, reproducible data and deeper scientific insights across disciplines.

Measure More



Capture 100+ phenotyping parameters – continuously, precisely, in real time

Trust the Data




Proven in 4000+ publications. PhenoMaster delivers reproducible data that stands up to any validation.

Stay in Sync



Every signal – metabolism, activity, behavior – synchronized in one continuous data stream.



400+

Trusted laboratories worldwide

4000+

Used in publications

15+

Integrated modules


Click to watch the most advanced metabolic phenotyping platform in action



From Metabolism to Behavior. Complete Phenotyping.

Benefit from the widest selection of modules available on the market – designed to deliver insights that no other system can provide.

ALL-IN-ONE PLATFORM INFINITE MODULES



Metabolism.

Reliable data. Reproducible results.



PhenoMaster NG captures every dimension of metabolic performance through continuous, synchronized data streams.

From oxygen consumption to respiratory exchange ratio, every parameter is measured with high temporal resolution and precision – allowing researchers to quantify energy balance, substrate utilization, and metabolic flexibility in real time.

The system integrates indirect calorimetry with advanced sensing of temperature, humidity, food and liquid intake, and body weight – all aligned in a unified data architecture for reliable and reproducible phenotyping.

ONE EXPERIMENT.
ONE TIMELINE.
COMPLETE SYNCHRONIZATION.

Long-term real-time monitoring of metabolic parameters in a cohort of mice shows precise synchronization across PhenoMaster modules.













1 Across alternating light and dark phases, PhenoMaster clearly captures the circadian rhythm of nocturnal rodents. Energy expenditure and activity patterns show perfect temporal alignment, reflecting real time synchronization

2 Drinking behavior follows the same rhythmic pattern – increasing at night alongside feeding and activity, and decreasing during the resting “light” phase.

3 During the fasting period (day 4), food intake stops and RER shifts and rebound accordingly, demonstrating metabolic adaptation captured by PhenoMaster.

CaloMax HS – The New Benchmark in Indirect Calorimetry

The most advanced indirect calorimetry technology ever engineered.

-  3× faster multiplex measurement engine
-  Dual-pump redundancy for continuous operation
-  Real-time hardware diagnostic for minimum downtime
-  Innovative Air-Conditioning system
-  Extended sensor lifespan for long-term reliability
-  Multi-sensor board for different gas sensors integration and streamline response
-  3 time smaller footprint
-  Automated moisture removal
-  Pre-filter unit for dust-free airflow
-  Minimum cabling for effortless handling

One platform.

Unlimited insight:




FLEXIBLE. SCALABLE. PRECISE.

PhenoMaster is a fully modular platform for indirect calorimetry and metabolic phenotyping – designed to evolve with your research.

Different cage architectures, gas analyzers, and sensor configurations can be flexibly combined to match any experimental design – from single-animal setups to large-scale high-throughput studies.

Continuous:	1 gas analyzer per cage, real-time 1s resolution	Short-term, acute studies
Multiplex:	Shared gas analyzers cycling between cages	Long-term, comparative experiments
Accelerated Multiplex:	Parallel analyzer pairs for simultaneous measurements	Cutting total runtime and expanding research flexibility

Beyond standard O₂ and CO₂ detection, the system also supports sensors for stable isotope tracing (¹³CO₂/¹²CO₂) and detection of CH₄ and H₂, expanding metabolic research into new analytical dimensions.



Behavior. Activity. Exercise.

PhenoMaster covers key aspects of rodent behavior – from feeding to cognition, from spontaneous activity to exercise.

The system enables both broad phenotyping and the ability to zoom in on specific behavioral features with unmatched precision.



INGESTIVE BEHAVIOR & BODY WEIGHT



- Precise sensors continuously track food, liquid intake and body weight
- Perfect for choice and preference studies
- Up to four access controls enable: scheduled, paired or yoked feeding, and deprivation paradigms

More details: 

Up to 4 Sensors per Cage

Spill & Leak Protection

High-Fat-Feeder

Programmable Access

RUNNING WHEELS



- Designed for long-term studies in standard, workload, and skilled versions
- Quantify exercise, motor skills, and circadian motor learning activity patterns
- Available in durable stainless steel or smooth closed-surface design

More details: 

Programmable Workload

Motor-Skilled Tests

Enable/Disable

Resistance Running

LOCOMOTOR ACTIVITY

ActiMot3



- Infrared activity frame with unparalleled 1.25 mm resolution
- From circadian patterns, to physical displacement and animals fine behavior
- Operates under any light conditions

Unique Cage Lid Holder

Finest Infrared Technology

5 mm Sensor Spacing

100+ Parameters

More details: 

OPERANT WALL

Cognition & Memory



- Unique module for cognitive and motivation testing in home cage
- Equipped with nose pokes, levers and sensors for long-term testing
- Ideal for extensive, long-term and fully-automated studies

Operant Behavior

Self-Administration


Food or Liquid Reward

Gut-Brain Axis

More details: 

WEIGHTLIFTING MODULE

Resistance Training



- Unique innovation for resistance training in rodents
- Adjustable weights enable standardized exercise paradigms
- Opens new perspectives in metabolism, aging, and muscle physiology research

Resistance Training

Muscle Adaptation

Adjustable Weights

Unique Innovation

More details: 

Unique Features

¹ MULTI-SENSOR GAS ANALYSIS

All gas sensors are integrated into one native multi-board from a single manufacturer – not a patchwork of devices. This ensures the fastest speed, the highest precision, and unmatched stability in metabolic readouts.

² SPECIAL CAGE LID HOLDER

Our exclusive design makes daily handling faster and easier. Bottles, feeders and bedding can be refilled in seconds – no need for extra hands or workspace. The built-in holder keeps the lid in place during animal handling, reducing workload and technical effort.

³ FAST CALIBRATION & DATA ANALYSIS



System calibration takes only 20 minutes before starting experiment. After data collection, results are processed quickly using the new Software Analytics Package – delivering comprehensive insights in record time.

⁴ ISOLATORS & ENVIRONMENTAL MONITORING



Keep your experiments clean and controlled. Cages can be placed inside isolators for germ-free or infectious studies, protecting both animals and researchers. The Lab Environmental Monitor continuously tracks CO₂, temperature, humidity, light and noise. Alarm system with push notifications ensures experiment safety.




Up to 18 cages per chamber



Transparent doors



Temperature, humidity, light control

5	CROSS-SPECIES COMBI SYSTEM	A single modular platform supporting both mice and rats. Adjustable airflow and shared core modules allow precise adaptation to different species and cage volumes – maximizing flexibility and efficiency.
6	GAS SENSOR FOR H₂ DETECTION	Optional gas sensor integrated into indirect calorimetry to measure gut-derived hydrogen production, enabling advanced analysis of microbiome activity.
7	EXTENSIVE CAGE SELECTION 	PhenoMaster offers a wide selection of specialized cages – including OptoCage, Urine & Feces Collection Cages, Basal Metabolic Rate Cage, Iso-Cage, Obese-Rat Cages, etc.
8	ADVANCED EXERCISE MODULES	Comprehensive suite of exercise tools – from voluntary running wheels to programmable resistance wheels and the Weightlifting Module.
9	CAGE-LEVEL CLIMATE MONITORING	Temperature and humidity sensors integrated into the cage lid track the exact microenvironment of each animal, ensuring consistent data and reproducible result.
10	HUMIDITY CONTROL	Humidity is a critical factor in metabolic research, influencing thermoregulation and energy expenditure. PhenoMaster is the only system allowing researchers to control this key variable inside the climate chamber.

SCIENTIFIC BREAKTHROUGHS POWERED BY PHENOMASTER



Discovering New Pathways in Metabolic Control

With the help of the PhenoMaster, Professor Matthias Tschöp and his team **identified a new class of multi-receptor drugs targeting GLP-1, GIP and glucagon pathways**, revolutionizing the treatment of obesity and metabolic diseases.



Pioneering Stable Isotope Integration


In 2018, TSE Systems pioneered the integration of ¹³C and ¹²C isotope sensors directly into the PhenoMaster indirect calorimetry system, enabling fully automated in vivo substrate tracing and **setting a new standard for metabolic research worldwide**.



Decoding the Microbiome in Controlled Environments

PhenoMaster platform enabled fully automated measurement of H₂ and CH₄ gases and introduced specialized isolators for germ-free animal studies, allowing researchers to investigate the complex dialogue between microbiota and host metabolism with unprecedented precision.




 Case study



 Articles



 Best practice


Analyze Faster. Discover More.

Turn complex datasets into clear insights with the new PhenoMaster Analytics Software, automated data preprocessing and standardized output significantly reduce analysis time – allowing you to focus on discovery instead of manual processing.




Powerful visualization tools and advanced statistics such as ANOVA, ANCOVA, and principal component analysis enable deep, flexible exploration of metabolic and behavioral results.

Built for reproducibility, collaboration, and speed – all in one seamless platform.

Investigation of energy expenditure and locomotor spontaneous activity in various animal cohort



Actogram Analysis




Name	Unit	De
CenA	[Cnts]	
CenF	[Cnts]	
CenT	[Cnts]	
CO2	[%]	
dCO2	[%]	
DistD	[cm]	
DistK	[cm]	
dO2	[%]	
Drink	[ml]	
DrinkC	[ml]	Drink (
Feed1	[g]	
Feed1C	[g]	Feed1 (
Flow	[l/min]	
H(1)	[kcal/h/kg]	
H(2)	[kcal/h/kg]	
H(3)	[kcal/h]	
Hum	[%]	
HumC	[%]	
LightC	[%]	
O2	[%]	
PerA	[Cnts]	
PerF	[Cnts]	
PerT	[Cnts]	
Press	[hPa]	
Ref.CO2		
Ref.O2		
Ref.SFlow		
RER		
S.Flow		
Speed		
Temp		
TempC		
VCO2(1)		




Assumptions check			
Univariate normality test			
Genotype	W	pval	normal
A	0.916	0.505	True
B	0.822	0.121	True
Homoscedasticity (equality of variance)			
	W	pval	equal_var
levene	0.056	0.82	True

ANCOVA and linear regression analysis of the relationship between EE and activity across genotypes



Statistical analysis across different genetic backgrounds and conditions



p2	p-tukey	eta-square
0.12	0.766	0.009
AN		

Why PhenoMaster?

All in One



MAXIMUM DATA FROM ONE EXPERIMENT

Over 15 integrated modules – including Calorimetry, Microbiome Gas Sensors, Activity, Feeding, Body Weight, Operant Wall, Running Wheel, Weight-Lifting Module, and Stable Isotope Tracing – deliver 100+ synchronized data readouts from a single setup.

Next Gen Technology



CONTINUOUS INNOVATION & SEAMLESS UPGRADES

Since 2004, PhenoMaster never stands still. Our team continuously develops new modules, sensors and analytics tools – keeping the platform at the forefront of research technology. Each upgrade integrates effortlessly, giving users early access to the latest innovation without replacing existing hardware.

Ethical Research



ENHANCED ANIMAL WELFARE & SPACE EFFICIENCY

Up to 30 % more engaged living space, reduced disturbance through non-invasive refilling, and a clean, open design that minimizes stress and intervention – because data quality starts with animal wellbeing.

Long Term Reliability




GERMAN QUALITY & SAFETY

Manufactured in Germany under strict ISO 9001 quality standards. Precision engineering ensures reliability, safety, and extended lifespan – benchmark of durability and trust. Dedicated technical & scientific support helps maintain performance and achieve consistent results.



The TSE PhenoMaster system has been critical to advancing our science, and our ability to precisely quantify a broad range of metabolic outcomes in our rodent models. The service team has aided us every step of the way in troubleshooting and keeping us up and running. Working with Gary has been amazing. He assisted us every step of the way from evaluating a purchase to ensuring prompt and attentive service even after the purchase.


Brian DeBosch
Associate Professor | Washington University
School of Medicine



Your next discovery starts here

- Connect

1 **CLICK HERE:**

A circular button with a teal background and a white mouse cursor icon pointing towards the bottom right.
- Consult

2 **BOOK A 20-MINUTE DISCOVERY CALL**

Discuss your research goals with our experts
- Configure

3 **GET YOUR TAILORED SYSTEM CONFIGURATION**

Optimized for your study design and parameters