

Product Overview

Sophisticated Life Science Research Instrumentation



TSE Dustless Precision Pellets

For research studies requiring precision feeding

www.TSE-Systems.com ■

Rel. March2005

◆ Dustless Precision Pellets

Features and advantages

- ◆ The pellets are nutritionally balanced and virtually dust-free.
- ◆ Hassle-free performance in automatic feeders.
- ◆ Suitable as a maintenance diet as well as a behavioral reward.
- ◆ Available in sizes of 14/ 20/ 45 mg.



- ◆ Tested routinely for hardness, weight and height.
- ◆ Produced in a GMP environment assuring accurate and consistent diet composition.
- ◆ For TSE pellet dispensers and other pellet dispensers.

Nutritional Assays

Nutritional assays and chemical screens are provided to ensure product consistency.

Ingredients

Protein, fat, fiber, ash and moisture

Purified Formula

Typical Nutritional Content

Protein	18.5%	Fat	4.9%
Fiber	4.5%	Ash	4.3%
Moisture	5.0%	Carbohydrate	61.5%

Caloric Value

3.62 cal/mg

Ingredients

Casein, sucrose, fiber, dextrose, DL-methionine, L-cystine, choline bitartrate, corn oil, corn syrup, mineral mix, vitamin mix, flowing and flavoring agents

Shelf Life

12 months in cool, dry storage

The nutritional content on every lot of diet is analyzed by an independent laboratory prior to shipment. The results of analysis are available upon request.

◆ Ordering Information

Cat-No.	Description
259901-PEL-14/50T	Precision-Pellets 14 mg, dustless BioServ Purified Formula For TSE Pellet Dispenser (14 mg) and other pellet dispensers. Pack = 50.000 pellets.
259901-PEL-20/50T	Precision-Pellets 20 mg, dustless BioServ Purified Formula For TSE Pellet Dispenser (20 mg) and other pellet dispensers. Pack = 50.000 pellets.
259901-PEL-45/50T	Precision-Pellets 45 mg, dustless BioServ Purified Formula For TSE Pellet Dispenser (45 mg) and other pellet dispensers. Pack = 50.000 pellets

■ Worldwide

TSE Systems GmbH
Siemensstrasse 21
61352 Bad Homburg
Germany

Phone: + 49-(0)6172-789-0
Fax: + 49-(0)6172-789-500

■ USA/Canada/Mexico

TSE Systems, Inc.
784 S. Poseyville Road
Midland, MI 48640
USA

Phone: 1-989-698-3067
Fax: 1-989-698-3068
Toll-Free (USA/Canada)
Phone: 1-866-466-8873
Fax: 1-866-467-8873

■ India

Axiom Biotek Inc.
Uniline House, 2nd Floor
198/23, Ramesh Market, East of Kailash
New Delhi 110 065
India

Phone: 0091-11-26469031
Fax: 0091-11-26481469
E-mail: harish@axiombiotek.com

www.TSE-Systems.com
info@TSE-Systems.com